



Myrian[®] XP-Prostate




Complete implementation of PI-RADS 2

- ◆ PI-RADS scoring mechanism integrated in the reading application
- ◆ Standard localization diagram implemented, providing valuable information for clinicians



Enhanced user experience

- ◆ Automatic sorting of sequences for accelerated loading
- ◆ Decoding of vendors private tags, allowing for exams handling from different manufacturers
- ◆ Multiple customizable hanging protocols for one-click operation, display of multiple sequences rearrangeable in a single click
- ◆ Compatible with all vendors or scanner generations 3 T, 1.5 T



Inspired
by patient
chain of care



Optimized communication with clinicians

- ◆ Unequaled structured reading, organ-based exhaustive reports
- ◆ Task-based workflow to report lesions, prostate & lymph nodes

intrasense.fr



Myrian[®]
XP-Prostate



intrasense

Dedicated reports

Customized
hanging protocol

Structured reading



Advanced
parametric maps



Multiparametric advanced visualization tools

- ◆ Specific time-saving tools available for ease-of-use
- ◆ Multiphase and multivolume tools: 3D & 2D display for T2, DCE and diffusion sequences
- ◆ Fusion of any sequence regardless of acquisition, orientation and thickness
- ◆ Built-in calculation of ADC, perfusion maps
- ◆ Time and intensity curve visualization built in the workflow



Broad clinical interest, usage, applications

- ◆ Diagnostics & scoring
- ◆ MRI guided biopsy
- ◆ Evaluation of response to treatment

CE 0459

FDA-
cleared